GuideLiner[®] Catheter

GuideLiner[®] Catheter Instructions For Use

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.


CAUTION

The GuideLiner catheter should be used by physicians with adequate training in the use of the device. The GuideLiner catheter is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infections. Contamination of the device may lead to illness or serious patient injury.

DEVICE DESCRIPTION

The GuideLiner is a single lumen catheter offered in three sizes compatible with 6F, 7F, and 8F guide catheters and may be placed over either an exchange length or 180cm guidewire. The 145cm device has a stainless steel shaft with a 20cm single lumen catheter of which the distal 15cm has a silicone wipe.

The GuideLiner catheter has a platinum-iridium marker band located 0.105" / 2.667mm from the distal tip, which enables visibility while using standard fluoroscopic methods. The device has two positioning marks located at 95cm (single mark) and 105cm (double mark) from the distal tip, respectively.

The GuideLiner catheters are delivered through a guide catheter resulting in an inner diameter that is approximately 1 French smaller than the guide catheter. The GuideLiner catheters are available in the following sizes and are compatible with guide catheters with the following sizes:

SPECIFICATIONS

Model	Compatible Guide Catheter	GuideLiner Min. I.D.	GuideLiner Tip O.D.
5571 6F (5-in-6)	≥ 6F (≥ 0.070" / 1.78mm I.D.)	0.056" / 1.42mm	0.067" / 1.70mm
5572 7F (6-in-7)	≥ 7F (≥ 0.078" / 1.98mm I.D.)	0.062" / 1.57mm	0.074" / 1.88mm
5573 8F (7-in-8)	≥ 8F (≥ 0.088" / 2.24mm I.D.)	0.071" / 1.80mm	0.084" / 2.13mm

INDICATIONS

GuideLiner catheters are intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement and exchange of guidewires and other interventional devices.

CONTRAINDICATIONS

The GuideLiner catheters are contraindicated in vessels less than 2.5mm in diameter, vessels in the neurovasculature and the venous system.

WARNINGS

Never advance the GuideLiner catheter into a vessel without a leading guidewire or without confirming location using fluoroscopic guidance. Vessel dissection or perforation may result.

Never advance the GuideLiner catheter into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the GuideLiner catheter, withdraw the GuideLiner catheter immediately.

Never advance the GuideLiner catheter more than 10cm beyond the tip of the guide catheter as the GuideLiner catheter may become lodged in the guide catheter making it difficult to remove.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

COMPLICATIONS

As with all catheterization procedures, complications with the GuideLiner catheter may occur. These may include:

- local or systemic infection
- air embolism
- intimal disruption
- arterial dissection
- perforation of the vessel wall
- vascular occlusion
- arterial thrombosis
- myocardial infarction
- arterial spasm

PRECAUTIONS

The GuideLiner catheter deployment procedure should be performed by physicians thoroughly trained in percutaneous intravascular techniques and procedures.

Do not use the GuideLiner catheter if the packaging has been damaged.

Inspect the GuideLiner catheter prior to use for any bends or kinks. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

The catheter lumen should be thoroughly flushed with heparinized saline prior to use to prevent clot formation.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solution should be considered.

Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the GuideLiner catheter. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

Each GuideLiner catheter includes:

- One GuideLiner catheter

Other materials required but not provided are:

- Guide catheter
- Y-adaptor with hemostasis valve (Tuohy-Borst type)
- Guidewire with diameter ≤ 0.014" / 0.36mm
- 10ml syringe for flushing
- Sterile heparinized saline (for system flushing)

PREPARATIONS FOR USE

1. Using sterile technique, transfer the dispenser coil with the GuideLiner catheter into the sterile field.
2. Flush the lumen of the GuideLiner catheter from the distal tip with heparinized saline solution.

DEPLOYMENT PROCEDURE

Deploy the GuideLiner catheter according to the following steps:

1. Secure the previously inserted guidewire and backload the distal tip of the GuideLiner catheter onto the guidewire and advance until it is just proximal to the hemostasis valve.
2. Open the hemostasis valve and advance the GuideLiner catheter through the hemostasis valve and into the guide catheter.
3. Under fluoroscopy, advance the GuideLiner catheter up to a maximum of 10cm beyond the distal end of the guide catheter and into the desired location within the vessel.

Warning: Never advance the GuideLiner catheter into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the GuideLiner catheter, withdraw the GuideLiner catheter immediately.

Warning: Never advance the GuideLiner catheter more than 10cm beyond the tip of the guide

English/Instructions for Use	1
Český/Návod k použití	2
Dansk/ Brugsanvisning	3
Nederlands/Gebruiksaanwijzing.....	4
Suomi/Käyttöohjeet.....	5
Français/Mode d'emploi	6
Deutsch/Gebrauchsanweisung	7
Ελληνικά/Οδηγίες χρήσης	8
Italiano/Istruzioni per l'uso.....	9
Norsk/Bruksanvisning	10
Polski/Instrukcja stosowania	11
Português/Instruções de Utilização.....	12
Español/Instrucciones de uso	13
Svenska/Bruksanvisning	14
Türkçe/Kullanım Talimatları	15

Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, MN 55369 USA
(888) 240-6001 USA
(763) 656-4300
(763) 656-4250
www.vasc.com

Vascular Solutions Zerus Limited
208 Business Innovation Centre
NUI Galway
Newcastle Road
Galway Ireland
+353 91 861611
+353 91 861612

catheter as the GuideLiner catheter may become lodged in the guide catheter making it difficult to remove.

- Using fluoroscopy, confirm the desired position of the GuideLiner catheter in the vessel.
- When inserted for an interventional procedure, position the metal collar of the GuideLiner catheter in a straight section of the guide catheter. Backload the interventional device over the in place guidewire and advance through the metal collar of the GuideLiner catheter into the desired vascular space.
- Tighten the Y-adaptor hemostasis valve securely on the proximal shaft of the GuideLiner catheter to prevent back-bleeding.
- Perform the catheterization procedure. After completing the procedure, remove the GuideLiner catheter prior to removing the guide catheter from the vessel.

PACKAGING & STORAGE

The GuideLiner catheter has been sterilized with ethylene oxide.

LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the GuideLiner catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special, or consequential damages arising from the use of the GuideLiner catheter. Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty.

No employee, agent, or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

International and U.S. patents pending.

GuideLiner® is a registered trademark of Vascular Solutions, Inc.

See the International Symbols Glossary on page 19.

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GuideLine® Vascular Solutions, Inc.'in tescilli bir ticari markasıdır.